


# MASTER\_MIX

350 nL in each well



Source

Destination


Plate Name: Source\_MATER\_MIX[1] << >>

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A																								
B																								
C																								
D																								
E																								
F																								
G																								
H																								
I																								
J																								
K																								
L																								
M																								
N																								
O																								
P																								


Plate Name: Destination\_MASTER\_MIX[1] << >>

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				
B	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				
C	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				
D	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				
E	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				
F	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				
G	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				
H	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				
I	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				
J	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1				
K	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
L	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
M	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
N	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
O	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
P	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1

[60  $\mu$ L in each well]



25, 50, 75 or 125 nL in different wells



Source

Destination

Plate Name: Source\_H2O[1] << >>

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A																								
B																								
C																								
D																								
E																								
F																								
G																								
H																								
I																								
J																								
K																								
L																								
M																								
N																								
O																								
P																								

Plate Name: Destination\_H2O[1] << >>

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A	1							1							1									
B	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
C	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
D	1							1							1									
E	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
F	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
G	1							1							1									
H	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
I	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
J	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
K	1							1							1							1	1	1
L	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
M	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
N	1							1							1							1	1	1
O	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
P	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1

**[60 µL in each well]**

**25 nL** (J2-J7, J9-J14, J18-J21, K22-K24, L22-L24, N22-N24, O22-O24, A1, A8, A15, D1, D8, D15, G1, G8, G15, K1, K8, K15, N1, N8, N15)

**50 nL** (B2-B7, B9-B14, B18-B21, E2-E7, E9-E14, E18-E21, H2-H7, H9-H14, H18-H21, L2-L7, L9-L14, L18-L21, O2-O7, O9-O14, O18-O21, M23, P23, J1, J8, J15)


**75 nL** (C2-C7, C9-C14, C18-C21, F2-F7, F9-F14, F18-F21, I2-I7, I9-I14, I18-I21, M2-M7, M9-M14, M18-M21, P2-P7, P9-P14, P18-P21, B1, B8, B15, E1, E8, E15, H1, H8, H15, L1, L8, L15, O1, O8, O15)

**100 nL** (C1, C8, C15, F1, F8, F15, I1, I8, I15, M1, M8, M15, P1, P8, P15)



# RNA

25 nL in each well



J1-J21, K22-24, L22-L24, M23,  
N22-N24, O22-O24, P23  
**= RNA- CONTROLS**

Source

Destination

Plate Name:

<<

>>

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	
A																									
B																									
C																									
D																									
E																									
F																									
G																									
H																									
I																									
J																									
K																									
L																									
M																									
N																									
O																									
P																									

**4  $\mu\text{g}/\mu\text{L}$  [L1] – 400 ng/ $\mu\text{L}$  [L2] – 40 ng/ $\mu\text{L}$  [L3]  
– 4 ng/ $\mu\text{L}$  [L4] – 400 pg/ $\mu\text{L}$  [L5]  
[20  $\mu\text{L}$  in each well]**

Plate Name:

<<


>>

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
B	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
C	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
D	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
E	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
F	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
G	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
H	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
I	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
J																								
K	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
L	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
M	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
N	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
O	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			
P	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1			

**100 ng [A1-I7] – 10 ng [A8-I14] – 1 ng [A16-  
I21] – 100 pg [K1-M21] – 10 pg [N1-P21]**

# RT\_PRIMERS

25 nL in each well



Source

Destination

Plate Name: Source\_RT\_PRIMERS[1] << >>

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A																								
B																								
C																								
D																								
E																								
F																								
G																								
H																								
I																								
J																								
K																								
L																								
M																								
N																								
O																								
P																								

**10  $\mu$ M [K1] – 20  $\mu$ M [K2] – 40  $\mu$ M [K3]  
– 80  $\mu$ M [K4] – 160  $\mu$ M [K5] – 240  $\mu$ M [K6]  
[20  $\mu$ L in each well]**

Plate Name: Destination\_RT\_PRIMERS[1] << >>

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
B		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
C		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
D		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
E		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
F		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
G		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
H		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
I		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
J		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
K		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
L		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
M		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
N		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
O		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				
P		1	1	1	1	1	1		1	1	1	1	1	1		1	1	1	1	1				

**0.5  $\mu$ M [col.2,9,16,K22,N22] – 1  $\mu$ M [col.  
3,10,17,K23,N23] – 2  $\mu$ M [col.4,11,18,K24, N24] –  
4  $\mu$ M [col.5,12,19,L22, O22] – 8  $\mu$ M [col.  
6,13,20,L23,O23] – 12  $\mu$ M [col. 7,14,21,L24, O24]**



# TSO (barcodes 1-70)


25, 50, or 100 nL (1, 2 or 4 drops)

A1-A7, D1-D7, G1-G7, J1-J21 → 4 drops  
 B1-B7, E1-E7, H1-H7 → 2 drops  
 C1-C7, F1-F7, I1-I7 → 1 drop

(K22-K24, L22-L24, M23 → 4 drops  
 N22-N24, O22-O24, P23 → 4 drops...)

Source

Destination



600  $\mu$ M (rows A,B,C)  
 400  $\mu$ M (rows D,E,F)  
 50  $\mu$ M (rows G,H,I,J)  
 [20  $\mu$ L in each well]

120  $\mu$ M (A1-A21, K1-K7, N1-N7) – 60  $\mu$ M (B1-B21, L1-L7, O1-O7) – 30  $\mu$ M (C1-C21, M1-M7, P1-P7) – 80  $\mu$ M (D1-D21, K8-K14, N8-N14, J1-J21, K22-K24, L22-L24, M23, N22-N24, O22-O24, P23) – 40  $\mu$ M (E1-E21, L8-L14, O8-O14) – 20  $\mu$ M (F1-F21, M8-M14, P8-P14) – 10  $\mu$ M (G1-G21, K15-K21, N15-N21) – 5  $\mu$ M (H1-H21, L15-L21, O15-O21) – 2.5  $\mu$ M (I1-I21, M15-M21, P1-P21)

# MASTER\_MIX PREPARATION (SSIII, nanoCAGE)

Reagent	Volume for 1 reaction (nL)	Stock conc.	Final conc.	Master_Mix for 384 reactions (201,25 µL)
Sorbitol/ Trehalose	40	0,66 M/3,3 M	0,0528M/0,264M	23
SuperScript III Reaction Buffer	100	5x	1x	57,5
DTT	50	0,1 M	0,01 M	28,75
dNTPs	31,25	10 mM	0,625 mM	17,97
Betain	75	5 M	0,75 M	43,13
SuperScript III	50	200 U/µL	20 U/uL	28,75
H <sub>2</sub> O	3,75	-	-	2,16
<b>TOTAL</b>	<b>350</b>			<b>201,25</b>

$(65 - 15) / 350 = 142$  destination wells filled per source well

→ 3x 65 µL wells required to effectively fill 384 wells

# MASTER\_MIX PREPARATION (SSIV, manufacturer)

Reagent	Volume for 1 reaction (nL)	Stock conc.	Final conc.	Master_Mix for 384 reactions (201,25 µL)
Sorbitol/ Trehalose	40	0,66 M/3,3 M	0,0528M/0,264M	23
SuperScript IV Reaction Buffer	100	5x	1x	57,5
DTT	25	0,1 M	0,005 M	14,375
dNTPs	31,25	10 mM	0,625 mM	17,97
Betain	75	5 M	0,75 M	43,13
SuperScript IV	25	200 U/µL	10 U/uL	14,375
H <sub>2</sub> O	53,75	-	-	30,90625
<b>TOTAL</b>	<b>350</b>			<b>201,25</b>

$(65 - 15) / 350 = 142$  destination wells filled per source well

→ 3x 65 µL wells required to effectively fill 384 wells

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24

A	1	2	3	4	5	6	7	1	2	3	4	5	6	7	1	2	3	4	5	6	7			
B	8	9	10	11	12	13	14	8	9	10	11	12	13	14	8	9	10	11	12	13	14			
C	15	16	17	18	19	20	21	15	16	17	18	19	20	21	15	16	17	18	19	20	21			
D	22	23	24	25	26	27	28	22	23	24	25	26	27	28	22	23	24	25	26	27	28			
E	29	30	31	32	33	34	35	29	30	31	32	33	34	35	29	30	31	32	33	34	35			
F	36	37	38	39	40	41	42	36	37	38	39	40	41	42	36	37	38	39	40	41	42			
G	43	44	45	46	47	48	49	43	44	45	46	47	48	49	43	44	45	46	47	48	49			
H	50	51	52	53	54	55	56	50	51	52	53	54	55	56	50	51	52	53	54	55	56			
I	57	58	59	60	61	62	63	57	58	59	60	61	62	63	57	58	59	60	61	62	63			
J	64	65	66	67	68	69	70	64	65	66	67	68	69	70	64	65	66	67	68	69	70			
K	1	2	3	4	5	6	7	22	23	24	25	26	27	28	43	44	45	46	47	48	49	65	66	67
L	8	9	10	11	12	13	14	29	30	31	32	33	34	35	50	51	52	53	54	55	56	68	69	70
M	15	16	17	18	19	20	21	36	37	38	39	40	41	42	57	58	59	60	61	62	63		64	
N	1	2	3	4	5	6	7	22	23	24	25	26	27	28	43	44	45	46	47	48	49	65	66	67
O	8	9	10	11	12	13	14	29	30	31	32	33	34	35	50	51	52	53	54	55	56	68	69	70
P	15	16	17	18	19	20	21	36	37	38	39	40	41	42	57	58	59	60	61	62	63		64	



		Random primer concentration (μM)						
		0	0,5	1	2	4	8	16
TSO  concentration  (μM)	160	no_RT_PRIMERS	320,000	160,000	80,000	40,000	20,000	10,000
	80	no_RT_PRIMERS	160,000	80,000	40,000	20,000	10,000	5,000
	40	no_RT_PRIMERS	80,000	40,000	20,000	10,000	5,000	2,500
	60	no_RT_PRIMERS	120,000	60,000	30,000	15,000	7,500	3,750
	30	no_RT_PRIMERS	60,000	30,000	15,000	7,500	3,750	1,875
	20	no_RT_PRIMERS	40,000	20,000	10,000	5,000	2,500	1,250
	10	no_RT_PRIMERS	20,000	10,000	5,000	2,500	1,250	0,625
	5	no_RT_PRIMERS	10,000	5,000	2,500	1,250	0,625	0,313
	2,5	no_RT_PRIMERS	5,000	2,500	1,250	0,625	0,313	0,156
	10	no_RT_PRIMERS	20,000	10,000	5,000	2,500	1,250	0,625

close to sc-nanoCAGE conditions

nanoCAGE conditions

no RNA

no RT primers

# RT REACTION

- ABI 7900 HT qPCR system used for the RT
- RT conditions SSIII (nanoCAGE):
  - 22°C, 10 min.
  - 50°C, 30 min.
  - 70°C, 15 min.
  - 4°C, Hold
- RT conditions SSIV (manufacturer):
  - 23°C, 10 min.
  - 50°C, 10 min.
  - 80°C, 10 min.
  - 4°C, Hold

# RT SAMPLES PURIFICATION

- 5  $\mu\text{L}$  of water added in 7 wells, then RT products from 70 wells were collected (5 RNA concentration tested)
- 2x AMPure XP purification : 1.2x
- Elutions in 30  $\mu\text{L}$  and 20  $\mu\text{L}$   $\text{H}_2\text{O}$



# qPCR

- Kapa Sybr Fast kit
- StepOne qPCR system
- Purified cDNA samples analysed in triplicates

# cDNA PCR

- Kapa HiFi Hot Start
- amplification cycles for each sample
- AMPure XP purification : 1x
- Elution in 20  $\mu\text{L}$   $\text{H}_2\text{O}$
- Picogreen quantitation and BioAnalyzer HS DNA chip