

Labcyte-RT Data Analysis (MOIRAI BED FILES)

Contents

| | |
|--|-----------|
| Data load and QC in R | 1 |
| Load CAGE libraries | 2 |
| Load summary statistics from MOIRAI and polish the names | 2 |
| Create a CAGEexp object and load expression data | 11 |
| Number of sequencing reads extracted per sample | 11 |
| CTSS ANALYSIS | 12 |
| Number of nanoCAGE tags mapping at CTSS positions in each group of samples | 12 |
| Number of nanoCAGE tags mapping at CTSS positions in each sample | 12 |
| Create CTSS clusters | 13 |
| Annotation with GENCODE | 13 |
| QC PLOTS | 14 |
| Boxplots | 14 |
| Dotplots | 25 |
| Counts | 27 |
| Barplots | 29 |
| Rarefaction | 30 |
| Correlation of CTSS clusters expression across samples (sample distance) | 32 |
| Create the correlation matrix | 32 |
| Heatmap | 32 |
| Principal Component Analysis (PCA) | 33 |
| Define PCA axis based on correlation matrix | 33 |
| PCA Plot 1: PC1 vs. PC2 annotated by PRIMERS_RATIO | 33 |
| Plate maps | 33 |
| RT primers | 33 |
| TSO | 34 |
| Ratio TSO / RT | 35 |
| Extracted | 36 |
| Artefacts | 37 |
| rDNA | 38 |
| Genes | 39 |
| Counts | 40 |
| Promoters | 41 |
| Grand summary plots | 42 |
| Extracted, artefacts, rDNA | 42 |
| Molecules, genes and promoters | 42 |

Data load and QC in R

```
BS_GENOME <- "BSgenome.Mmusculus.UCSC.mm9"
library(BS_GENOME, character.only = T)
```

```
library(CAGEr)
library(data.table)
library(ggplot2)
library(gplots)
library('RColorBrewer')
library(magrittr)
library(plyr)
library(MultiAssayExperiment)
library(SummarizedExperiment)
library(reshape)
library(vegan)
```

MOIRAI shortcuts

```
WORKFLOW      <- "OP-WORKFLOW-CAGEscan-short-reads-v2.1-rc1"
MOIRAI_PROJ    <- "project/Labcyte"
MOIRAI_USER    <- "nanoCAGE2"
ASSEMBLY      <- "mm9"
BASEDIR       <- "/osc-fs_home/scratch/moirai"
MOIRAI_BASE    <- file.path(BASEDIR, MOIRAI_USER)
MOIRAI_RESULTS <- file.path(MOIRAI_BASE, MOIRAI_PROJ)
```

Load CAGE libraries

Load summary statistics from MOIRAI and polish the names

```
libs <- smallCAGEqc::loadMoiraiStats(
  pipeline = "OP-WORKFLOW-CAGEscan-short-reads-v2.0",
  multiplex = file.path(MOIRAI_BASE, "input/171227_M00528_0321_000000000-B4GLP.multiplex.txt"),
  summary = file.path(MOIRAI_RESULTS, "171227_M00528_0321_000000000-B4GLP.OP-WORKFLOW-CAGEscan-short-")
libs$barcode_ID <- c(1:70)
libs$inputFiles <- list.files(path = "/osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP",
  pattern = "*.bed", full.names = TRUE)
libs$inputFiles <- paste0("/osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP",
  libs$inputFiles, ".bed")
libs$inputFileType <- c("bed")
libs$sampleLabels <- as.character(libs$samplename)
#libs

plate <- read.table("plate.txt", sep = "\t", header = TRUE)
plate_10ng <- subset(plate, plate$RNA == 10)
plate_10ng_no_RNA <- plate[224:230,]
plate_10ng_all <- rbind(plate_10ng, plate_10ng_no_RNA)
#plate_10ng_all

libs <- cbind(libs, plate_10ng_all)
libs[,24] <- NULL
#rownames(libs) <- NULL
libs$PRIMERS_RATIO <- sub("no_RT_PRIMERS", "NA", libs$PRIMERS_RATIO)
libs$PRIMERS_RATIO <- as.numeric(libs$PRIMERS_RATIO)

## Warning: NAs introduced by coercion

libs$RNA <- as.numeric(libs$RNA)
libs <- libs[order(-libs$PRIMERS_RATIO),]
```

```

libs <- libs[order(~libs$RNA),]
libs$PRIMERS_RATIO[is.na(libs$PRIMERS_RATIO)] <- "no_RT_PRIMERS"
libs

```

| ## | samplename | group | barcode | index | total | extracted | cleaned | tagdust |
|----|------------|--------|---------|----------|-------|-----------|---------|---------|
| ## | ACACGT | ACACGT | ACACGT | TAGGCATG | 0 | 39556 | 30362 | 3435 |
| ## | AACTC | AACTC | AACTC | TAGGCATG | 0 | 159844 | 121202 | 16681 |
| ## | ACATGA | ACATGA | ACATGA | TAGGCATG | 0 | 32474 | 25537 | 3094 |
| ## | ACAGAT | ACAGAT | ACAGAT | TAGGCATG | 0 | 138710 | 109459 | 14604 |
| ## | AGTACG | AGTACG | AGTACG | TAGGCATG | 0 | 14180 | 10760 | 1275 |
| ## | AGTGAT | AGTGAT | AGTGAT | TAGGCATG | 0 | 168934 | 137231 | 12654 |
| ## | ACAGCA | ACAGCA | ACAGCA | TAGGCATG | 0 | 158662 | 118015 | 17660 |
| ## | AGTAGC | AGTAGC | AGTAGC | TAGGCATG | 0 | 56598 | 43154 | 4590 |
| ## | AGTGCA | AGTGCA | AGTGCA | TAGGCATG | 0 | 428768 | 340280 | 36348 |
| ## | ATCGCA | ATCGCA | ATCGCA | TAGGCATG | 0 | 43925 | 30747 | 4033 |
| ## | ACAGTG | ACAGTG | ACAGTG | TAGGCATG | 0 | 218169 | 162650 | 24469 |
| ## | AGTATA | AGTATA | AGTATA | TAGGCATG | 0 | 15195 | 10726 | 1692 |
| ## | AGTGTG | AGTGTG | AGTGTG | TAGGCATG | 0 | 16988 | 11880 | 1582 |
| ## | ATCGTG | ATCGTG | ATCGTG | TAGGCATG | 0 | 5852 | 3888 | 545 |
| ## | CACATA | CACATA | CACATA | TAGGCATG | 0 | 49911 | 31480 | 6737 |
| ## | ACATAC | ACATAC | ACATAC | TAGGCATG | 0 | 252993 | 170605 | 43474 |
| ## | AGTCAG | AGTCAG | AGTCAG | TAGGCATG | 0 | 18231 | 12169 | 2451 |
| ## | ATCACG | ATCACG | ATCACG | TAGGCATG | 0 | 11697 | 7038 | 1808 |
| ## | ATCTAC | ATCTAC | ATCTAC | TAGGCATG | 0 | 61931 | 39776 | 7975 |
| ## | CACGAT | CACGAT | CACGAT | TAGGCATG | 0 | 14078 | 8046 | 1768 |
| ## | CGACAG | CGACAG | CGACAG | TAGGCATG | 0 | 5787 | 3355 | 763 |
| ## | AGTCGT | AGTCGT | AGTCGT | TAGGCATG | 0 | 42894 | 27736 | 5143 |
| ## | ATCAGC | ATCAGC | ATCAGC | TAGGCATG | 0 | 1053249 | 719844 | 99891 |
| ## | ATCTCT | ATCTCT | ATCTCT | TAGGCATG | 0 | 7558 | 4713 | 951 |
| ## | CACGCA | CACGCA | CACGCA | TAGGCATG | 0 | 168925 | 101801 | 20625 |
| ## | CGACGT | CGACGT | CGACGT | TAGGCATG | 0 | 50902 | 28765 | 5996 |
| ## | CGATCT | CGATCT | CGATCT | TAGGCATG | 0 | 26817 | 16094 | 3128 |
| ## | ATCATA | ATCATA | ATCATA | TAGGCATG | 0 | 17945 | 10155 | 3332 |
| ## | ATCTGA | ATCTGA | ATCTGA | TAGGCATG | 0 | 198657 | 131168 | 17902 |
| ## | CACGTG | CACGTG | CACGTG | TAGGCATG | 0 | 30261 | 17471 | 3498 |
| ## | CGACTC | CGACTC | CGACTC | TAGGCATG | 0 | 151581 | 86422 | 16624 |
| ## | CGATGA | CGATGA | CGATGA | TAGGCATG | 0 | 3492 | 2084 | 483 |
| ## | CTGCTC | CTGCTC | CTGCTC | TAGGCATG | 0 | 12831 | 7704 | 1544 |
| ## | CACACG | CACACG | CACACG | TAGGCATG | 0 | 1333918 | 777491 | 188522 |
| ## | CACTAC | CACTAC | CACTAC | TAGGCATG | 0 | 551765 | 297147 | 87676 |
| ## | CGAGAT | CGAGAT | CGAGAT | TAGGCATG | 0 | 270601 | 151394 | 29998 |
| ## | CTGACG | CTGACG | CTGACG | TAGGCATG | 0 | 15719 | 8510 | 2101 |
| ## | CTGTAC | CTGTAC | CTGTAC | TAGGCATG | 0 | 7689 | 4464 | 1251 |
| ## | GAGCAG | GAGCAG | GAGCAG | TAGGCATG | 0 | 11251 | 6691 | 1313 |
| ## | CACTCT | CACTCT | CACTCT | TAGGCATG | 0 | 399697 | 210815 | 58589 |
| ## | CGAGCA | CGAGCA | CGAGCA | TAGGCATG | 0 | 135396 | 69981 | 18350 |
| ## | CTGAGC | CTGAGC | CTGAGC | TAGGCATG | 0 | 31942 | 17171 | 4917 |
| ## | CTGTCT | CTGTCT | CTGTCT | TAGGCATG | 0 | 231518 | 60117 | 145048 |
| ## | GAGCGT | GAGCGT | GAGCGT | TAGGCATG | 0 | 11654 | 1146 | 245 |
| ## | CGAGTG | CGAGTG | CGAGTG | TAGGCATG | 0 | 43603 | 20443 | 6994 |
| ## | CTGATA | CTGATA | CTGATA | TAGGCATG | 0 | 45389 | 22937 | 8017 |
| ## | CTGTGA | CTGTGA | CTGTGA | TAGGCATG | 0 | 29711 | 15875 | 4364 |
| ## | GAGCTC | GAGCTC | GAGCTC | TAGGCATG | 0 | 19222 | 11463 | 1980 |
| ## | CTGCAG | CTGCAG | CTGCAG | TAGGCATG | 0 | 92305 | 47128 | 13132 |

| | | | | | | | | | |
|----|--------|--------|--------|--------|-------------|--------|------------|-------|-------|
| ## | GAGACG | GAGACG | GAGACG | GAGACG | TAGGCATG | 0 | 10894 | 5736 | 1614 |
| ## | GAGTAC | GAGTAC | GAGTAC | GAGTAC | TAGGCATG | 0 | 917 | 492 | 227 |
| ## | GAGAGC | GAGAGC | GAGAGC | GAGAGC | TAGGCATG | 0 | 92220 | 44459 | 11116 |
| ## | GAGTCT | GAGTCT | GAGTCT | GAGTCT | TAGGCATG | 0 | 2006 | 1237 | 312 |
| ## | GAGTGA | GAGTGA | GAGTGA | GAGTGA | TAGGCATG | 0 | 17132 | 8698 | 2693 |
| ## | ACACAG | ACACAG | ACACAG | ACACAG | TAGGCATG | 0 | 42054 | 32450 | 4446 |
| ## | ACATCT | ACATCT | ACATCT | ACATCT | TAGGCATG | 0 | 33745 | 27224 | 3571 |
| ## | AGTCTC | AGTCTC | AGTCTC | AGTCTC | TAGGCATG | 0 | 3059 | 2204 | 193 |
| ## | ATCGAT | ATCGAT | ATCGAT | ATCGAT | TAGGCATG | 0 | 10403 | 7116 | 762 |
| ## | CACAGC | CACAGC | CACAGC | CACAGC | TAGGCATG | 0 | 8329 | 6011 | 604 |
| ## | CACTGA | CACTGA | CACTGA | CACTGA | TAGGCATG | 0 | 6506 | 3906 | 712 |
| ## | CGATAC | CGATAC | CGATAC | CGATAC | TAGGCATG | 0 | 1633 | 1032 | 117 |
| ## | CTGCGT | CTGCGT | CTGCGT | CTGCGT | TAGGCATG | 0 | 1318 | 973 | 67 |
| ## | GAGATA | GAGATA | GAGATA | GAGATA | TAGGCATG | 0 | 717 | 435 | 54 |
| ## | GCTAGC | GCTAGC | GCTAGC | GCTAGC | TAGGCATG | 0 | 1895 | 1202 | 279 |
| ## | GCTATA | GCTATA | GCTATA | GCTATA | TAGGCATG | 0 | 1206 | 786 | 219 |
| ## | GCTCAG | GCTCAG | GCTCAG | GCTCAG | TAGGCATG | 0 | 2855 | 1629 | 829 |
| ## | GCTCGT | GCTCGT | GCTCGT | GCTCGT | TAGGCATG | 0 | 3098 | 1593 | 1061 |
| ## | GCTCTC | GCTCTC | GCTCTC | GCTCTC | TAGGCATG | 0 | 4055 | 1809 | 1996 |
| ## | GCTGAT | GCTGAT | GCTGAT | GCTGAT | TAGGCATG | 0 | 23130 | 8235 | 13303 |
| ## | GCTACG | GCTACG | GCTACG | GCTACG | TAGGCATG | 0 | 915 | 584 | 71 |
| ## | | rdna | spikes | mapped | properpairs | counts | barcode_ID | | |
| ## | ACACGT | 5752 | 7 | 24315 | 19174 | 7757 | 2 | | |
| ## | ACACTC | 21910 | 51 | 97291 | 79126 | 33749 | 3 | | |
| ## | ACATGA | 3834 | 9 | 19931 | 16201 | 6915 | 9 | | |
| ## | ACAGAT | 14608 | 39 | 86021 | 70039 | 23945 | 4 | | |
| ## | AGTACG | 2142 | 3 | 7677 | 6221 | 2515 | 10 | | |
| ## | AGTGAT | 18983 | 66 | 96297 | 79334 | 39870 | 16 | | |
| ## | ACAGCA | 22924 | 63 | 93817 | 76382 | 20614 | 5 | | |
| ## | AGTAGC | 8832 | 22 | 30252 | 24577 | 8761 | 11 | | |
| ## | AGTGCA | 51993 | 146 | 247202 | 203106 | 73855 | 17 | | |
| ## | ATCGCA | 9134 | 11 | 22802 | 18594 | 6957 | 23 | | |
| ## | ACAGTG | 30957 | 93 | 126928 | 103827 | 30539 | 6 | | |
| ## | AGTATA | 2763 | 14 | 7293 | 5944 | 2650 | 12 | | |
| ## | AGTGTG | 3523 | 3 | 8006 | 6661 | 2843 | 18 | | |
| ## | ATCGTG | 1417 | 2 | 2848 | 2226 | 881 | 24 | | |
| ## | CACATA | 11680 | 14 | 25125 | 20627 | 6761 | 30 | | |
| ## | ACATAC | 38832 | 82 | 131967 | 109128 | 28385 | 7 | | |
| ## | AGTCAG | 3607 | 4 | 8717 | 7070 | 2536 | 13 | | |
| ## | ATCACG | 2847 | 4 | 4928 | 3893 | 1880 | 19 | | |
| ## | ATCTAC | 14171 | 9 | 27992 | 22593 | 8438 | 25 | | |
| ## | CACGAT | 4261 | 3 | 6151 | 5052 | 1738 | 31 | | |
| ## | CGACAG | 1665 | 4 | 2642 | 2102 | 614 | 37 | | |
| ## | AGTCGT | 10008 | 7 | 19690 | 16325 | 4776 | 14 | | |
| ## | ATCAGC | 233207 | 306 | 526082 | 438141 | 138913 | 20 | | |
| ## | ATCTCT | 1892 | 2 | 3025 | 2277 | 954 | 26 | | |
| ## | CACGCA | 46442 | 57 | 80641 | 65816 | 12532 | 32 | | |
| ## | CGACGT | 16127 | 14 | 22543 | 18241 | 4303 | 38 | | |
| ## | CGATCT | 7592 | 3 | 12355 | 9934 | 2887 | 44 | | |
| ## | ATCATA | 4454 | 4 | 6651 | 5320 | 2136 | 21 | | |
| ## | ATCTGA | 49511 | 76 | 90423 | 74762 | 24103 | 27 | | |
| ## | CACGTG | 9282 | 10 | 13699 | 11134 | 2336 | 33 | | |
| ## | CGACTC | 48495 | 40 | 67288 | 55103 | 9963 | 39 | | |
| ## | CGATGA | 924 | 1 | 1537 | 1181 | 368 | 45 | | |

| | | | | | | | |
|----|--------|---|-----|--------|--------|-------|----|
| ## | CTGCTC | 3583 | 0 | 5979 | 4164 | 1424 | 51 |
| ## | CACACG | 367550 | 355 | 615543 | 505838 | 97382 | 28 |
| ## | CACTAC | 166824 | 118 | 230182 | 188588 | 37193 | 34 |
| ## | CGAGAT | 89139 | 70 | 116804 | 95921 | 19089 | 40 |
| ## | CTGACG | 5100 | 8 | 6343 | 5141 | 1441 | 46 |
| ## | CTGTAC | 1972 | 2 | 3451 | 2486 | 651 | 52 |
| ## | GAGCAG | 3239 | 8 | 4988 | 3427 | 1223 | 58 |
| ## | CACTCT | 130193 | 100 | 162843 | 132950 | 27316 | 35 |
| ## | CGAGCA | 47019 | 45 | 54182 | 44151 | 8092 | 41 |
| ## | CTGAGC | 9841 | 12 | 13201 | 10507 | 2499 | 47 |
| ## | CTGTCT | 26325 | 28 | 42593 | 30380 | 7757 | 53 |
| ## | GAGCGT | 10262 | 1 | 787 | 519 | 202 | 59 |
| ## | CGAGTG | 16144 | 10 | 15140 | 12073 | 2785 | 42 |
| ## | CTGATA | 14419 | 16 | 16897 | 13603 | 3579 | 48 |
| ## | CTGTGA | 9462 | 10 | 11972 | 9535 | 2790 | 54 |
| ## | GAGCTC | 5776 | 3 | 8345 | 6477 | 2030 | 60 |
| ## | CTGCAG | 32013 | 32 | 36642 | 29777 | 5820 | 49 |
| ## | GAGACG | 3543 | 1 | 3956 | 3074 | 1121 | 55 |
| ## | GAGTAC | 197 | 1 | 316 | 235 | 106 | 61 |
| ## | GAGAGC | 36627 | 18 | 30852 | 25105 | 7384 | 56 |
| ## | GAGTCT | 457 | 0 | 786 | 575 | 231 | 62 |
| ## | GAGTGA | 5738 | 3 | 5917 | 4594 | 1647 | 63 |
| ## | ACACAG | 5146 | 12 | 26243 | 20877 | 11074 | 1 |
| ## | ACATCT | 2945 | 5 | 21474 | 16282 | 8834 | 8 |
| ## | AGTCTC | 657 | 5 | 1542 | 1200 | 756 | 15 |
| ## | ATCGAT | 2523 | 2 | 5221 | 4391 | 3043 | 22 |
| ## | CACAGC | 1710 | 4 | 4673 | 3616 | 1346 | 29 |
| ## | CACTGA | 1887 | 1 | 3069 | 2420 | 1370 | 36 |
| ## | CGATAC | 483 | 1 | 581 | 339 | 166 | 43 |
| ## | CTGCGT | 278 | 0 | 640 | 396 | 170 | 50 |
| ## | GAGATA | 228 | 0 | 276 | 200 | 108 | 57 |
| ## | GCTAGC | 414 | 0 | 815 | 657 | 454 | 65 |
| ## | GCTATA | 201 | 0 | 512 | 371 | 206 | 66 |
| ## | GCTCAG | 397 | 0 | 1154 | 848 | 536 | 67 |
| ## | GCTCGT | 442 | 2 | 1001 | 798 | 546 | 68 |
| ## | GCTCTC | 249 | 1 | 1214 | 571 | 231 | 69 |
| ## | GCTGAT | 1581 | 11 | 3721 | 2785 | 2054 | 70 |
| ## | GCTACG | 260 | 0 | 437 | 367 | 281 | 64 |
| ## | | | | | | | |
| ## | ACACGT | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | ACACTC | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | ACATGA | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | ACAGAT | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | AGTACG | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | AGTGAT | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | ACAGCA | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | AGTAGC | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | AGTGCA | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | ATCGCA | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | ACAGTG | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | AGTATA | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | AGTGTG | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | ATCGTG | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |
| ## | CACATA | /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W | | | | | |


```

## GCTACG /osc-fs_home/scratch/moirai/nanoCAGE2/project/Labcyte/171227_M00528_0321_000000000-B4GLP.OP-W
##      inputFileType sampleLabels Well Row Col MASTER_MIX_vol      TSD
## ACACGT          bed      ACACGT  A09  A   9              350 80.0000
## AACTC           bed      AACTC   A10  A  10              350 80.0000
## ACATGA          bed      ACATGA  B09  B   9              350 40.0000
## ACAGAT          bed      ACAGAT  A11  A  11              350 80.0000
## AGTACG          bed      AGTACG  B10  B  10              350 40.0000
## AGTGAT          bed      AGTGAT  C09  C   9              350 20.0000
## ACAGCA          bed      ACAGCA  A12  A  12              350 80.0000
## AGTAGC          bed      AGTAGC  B11  B  11              350 40.0000
## AGTGCA          bed      AGTGCA  C10  C  10              350 20.0000
## ATCGCA          bed      ATCGCA  D09  D   9              350 10.0000
## ACAGTG          bed      ACAGTG  A13  A  13              350 80.0000
## AGTATA          bed      AGTATA  B12  B  12              350 40.0000
## AGTGTG          bed      AGTGTG  C11  C  11              350 20.0000
## ATCGTG          bed      ATCGTG  D10  D  10              350 10.0000
## CACATA          bed      CACATA  E09  E   9              350  5.0000
## ACATAC          bed      ACATAC  A14  A  14              350 80.0000
## AGTCAG          bed      AGTCAG  B13  B  13              350 40.0000
## ATCACG          bed      ATCACG  C12  C  12              350 20.0000
## ATCTAC          bed      ATCTAC  D11  D  11              350 10.0000
## CACGAT          bed      CACGAT  E10  E  10              350  5.0000
## CGACAG          bed      CGACAG  F09  F   9              350  2.5000
## AGTCGT          bed      AGTCGT  B14  B  14              350 40.0000
## ATCAGC          bed      ATCAGC  C13  C  13              350 20.0000
## ATCTCT          bed      ATCTCT  D12  D  12              350 10.0000
## CACGCA          bed      CACGCA  E11  E  11              350  5.0000
## CGACGT          bed      CGACGT  F10  F  10              350  2.5000
## CGATCT          bed      CGATCT  G09  G   9              350  1.2500
## ATCATA          bed      ATCATA  C14  C  14              350 20.0000
## ATCTGA          bed      ATCTGA  D13  D  13              350 10.0000
## CACGTG          bed      CACGTG  E12  E  12              350  5.0000
## CGACTC          bed      CGACTC  F11  F  11              350  2.5000
## CGATGA          bed      CGATGA  G10  G  10              350  1.2500
## CTGCTC          bed      CTGCTC  H09  H   9              350  0.6250
## CACACG          bed      CACACG  D14  D  14              350 10.0000
## CACTAC          bed      CACTAC  E13  E  13              350  5.0000
## CGAGAT          bed      CGAGAT  F12  F  12              350  2.5000
## CTGACG          bed      CTGACG  G11  G  11              350  1.2500
## CTGTAC          bed      CTGTAC  H10  H  10              350  0.6250
## GAGCAG          bed      GAGCAG  I09  I   9              350  0.3125
## CACTCT          bed      CACTCT  E14  E  14              350  5.0000
## CGAGCA          bed      CGAGCA  F13  F  13              350  2.5000
## CTGAGC          bed      CTGAGC  G12  G  12              350  1.2500
## CTGTCT          bed      CTGTCT  H11  H  11              350  0.6250
## GAGCGT          bed      GAGCGT  I10  I  10              350  0.3125
## CGAGTG          bed      CGAGTG  F14  F  14              350  2.5000
## CTGATA          bed      CTGATA  G13  G  13              350  1.2500
## CTGTGA          bed      CTGTGA  H12  H  12              350  0.6250
## GAGCTC          bed      GAGCTC  I11  I  11              350  0.3125
## CTGCAG          bed      CTGCAG  G14  G  14              350  1.2500
## GAGACG          bed      GAGACG  H13  H  13              350  0.6250
## GAGTAC          bed      GAGTAC  I12  I  12              350  0.3125
## GAGAGC          bed      GAGAGC  H14  H  14              350  0.6250

```

| | | | | | | | | |
|----|---------|------------|----------------|-----|---------|---------|--------------|---------|
| ## | GAGTCT | bed | GAGTCT | I13 | I | 13 | 350 | 0.3125 |
| ## | GAGTGA | bed | GAGTGA | I14 | I | 14 | 350 | 0.3125 |
| ## | ACACAG | bed | ACACAG | A08 | A | 8 | 350 | 80.0000 |
| ## | ACATCT | bed | ACATCT | B08 | B | 8 | 350 | 40.0000 |
| ## | AGTCTC | bed | AGTCTC | C08 | C | 8 | 350 | 20.0000 |
| ## | ATCGAT | bed | ATCGAT | D08 | D | 8 | 350 | 10.0000 |
| ## | CACAGC | bed | CACAGC | E08 | E | 8 | 350 | 5.0000 |
| ## | CACTGA | bed | CACTGA | F08 | F | 8 | 350 | 2.5000 |
| ## | CGATAC | bed | CGATAC | G08 | G | 8 | 350 | 1.2500 |
| ## | CTGCGT | bed | CTGCGT | H08 | H | 8 | 350 | 0.6250 |
| ## | GAGATA | bed | GAGATA | I08 | I | 8 | 350 | 0.3125 |
| ## | GCTAGC | bed | GCTAGC | J09 | J | 9 | 350 | 10.0000 |
| ## | GCTATA | bed | GCTATA | J10 | J | 10 | 350 | 10.0000 |
| ## | GCTCAG | bed | GCTCAG | J11 | J | 11 | 350 | 10.0000 |
| ## | GCTCGT | bed | GCTCGT | J12 | J | 12 | 350 | 10.0000 |
| ## | GCTCTC | bed | GCTCTC | J13 | J | 13 | 350 | 10.0000 |
| ## | GCTGAT | bed | GCTGAT | J14 | J | 14 | 350 | 10.0000 |
| ## | GCTACG | bed | GCTACG | J08 | J | 8 | 350 | 10.0000 |
| ## | TSO_vol | RT_PRIMERS | RT_PRIMERS_vol | RNA | RNA_vol | H2O_vol | total_volume | |
| ## | ACACGT | 100 | 0.125 | 25 | 10 | 25 | 0 | 500 |
| ## | AACTC | 100 | 0.250 | 25 | 10 | 25 | 0 | 500 |
| ## | ACATGA | 50 | 0.125 | 25 | 10 | 25 | 50 | 500 |
| ## | ACAGAT | 100 | 0.500 | 25 | 10 | 25 | 0 | 500 |
| ## | AGTACG | 50 | 0.250 | 25 | 10 | 25 | 50 | 500 |
| ## | AGTGAT | 25 | 0.125 | 25 | 10 | 25 | 75 | 500 |
| ## | ACAGCA | 100 | 1.000 | 25 | 10 | 25 | 0 | 500 |
| ## | AGTAGC | 50 | 0.500 | 25 | 10 | 25 | 50 | 500 |
| ## | AGTGCA | 25 | 0.250 | 25 | 10 | 25 | 75 | 500 |
| ## | ATCGCA | 100 | 0.125 | 25 | 10 | 25 | 0 | 500 |
| ## | ACAGTG | 100 | 2.000 | 25 | 10 | 25 | 0 | 500 |
| ## | AGTATA | 50 | 1.000 | 25 | 10 | 25 | 50 | 500 |
| ## | AGTGTG | 25 | 0.500 | 25 | 10 | 25 | 75 | 500 |
| ## | ATCGTG | 100 | 0.250 | 25 | 10 | 25 | 0 | 500 |
| ## | CACATA | 50 | 0.125 | 25 | 10 | 25 | 50 | 500 |
| ## | ACATAC | 100 | 4.000 | 25 | 10 | 25 | 0 | 500 |
| ## | AGTCAG | 50 | 2.000 | 25 | 10 | 25 | 50 | 500 |
| ## | ATCACG | 25 | 1.000 | 25 | 10 | 25 | 75 | 500 |
| ## | ATCTAC | 100 | 0.500 | 25 | 10 | 25 | 0 | 500 |
| ## | CACGAT | 50 | 0.250 | 25 | 10 | 25 | 50 | 500 |
| ## | CGACAG | 25 | 0.125 | 25 | 10 | 25 | 75 | 500 |
| ## | AGTCGT | 50 | 4.000 | 25 | 10 | 25 | 50 | 500 |
| ## | ATCAGC | 25 | 2.000 | 25 | 10 | 25 | 75 | 500 |
| ## | ATCTCT | 100 | 1.000 | 25 | 10 | 25 | 0 | 500 |
| ## | CACGCA | 50 | 0.500 | 25 | 10 | 25 | 50 | 500 |
| ## | CGACGT | 25 | 0.250 | 25 | 10 | 25 | 75 | 500 |
| ## | CGATCT | 100 | 0.125 | 25 | 10 | 25 | 0 | 500 |
| ## | ATCATA | 25 | 4.000 | 25 | 10 | 25 | 75 | 500 |
| ## | ATCTGA | 100 | 2.000 | 25 | 10 | 25 | 0 | 500 |
| ## | CACGTG | 50 | 1.000 | 25 | 10 | 25 | 50 | 500 |
| ## | CGACTC | 25 | 0.500 | 25 | 10 | 25 | 75 | 500 |
| ## | CGATGA | 100 | 0.250 | 25 | 10 | 25 | 0 | 500 |
| ## | CTGCTC | 50 | 0.125 | 25 | 10 | 25 | 50 | 500 |
| ## | CACACG | 100 | 4.000 | 25 | 10 | 25 | 0 | 500 |
| ## | CACTAC | 50 | 2.000 | 25 | 10 | 25 | 50 | 500 |

| | | | | | | | | |
|----|---------------|-----|-------|----|----|----|-----|-----|
| ## | CGAGAT | 25 | 1.000 | 25 | 10 | 25 | 75 | 500 |
| ## | CTGACG | 100 | 0.500 | 25 | 10 | 25 | 0 | 500 |
| ## | CTGTAC | 50 | 0.250 | 25 | 10 | 25 | 50 | 500 |
| ## | GAGCAG | 25 | 0.125 | 25 | 10 | 25 | 75 | 500 |
| ## | CACTCT | 50 | 4.000 | 25 | 10 | 25 | 50 | 500 |
| ## | CGAGCA | 25 | 2.000 | 25 | 10 | 25 | 75 | 500 |
| ## | CTGAGC | 100 | 1.000 | 25 | 10 | 25 | 0 | 500 |
| ## | CTGTCT | 50 | 0.500 | 25 | 10 | 25 | 50 | 500 |
| ## | GAGCGT | 25 | 0.250 | 25 | 10 | 25 | 75 | 500 |
| ## | CGAGTG | 25 | 4.000 | 25 | 10 | 25 | 75 | 500 |
| ## | CTGATA | 100 | 2.000 | 25 | 10 | 25 | 0 | 500 |
| ## | CTGTGA | 50 | 1.000 | 25 | 10 | 25 | 50 | 500 |
| ## | GAGCTC | 25 | 0.500 | 25 | 10 | 25 | 75 | 500 |
| ## | CTGCAG | 100 | 4.000 | 25 | 10 | 25 | 0 | 500 |
| ## | GAGACG | 50 | 2.000 | 25 | 10 | 25 | 50 | 500 |
| ## | GAGTAC | 25 | 1.000 | 25 | 10 | 25 | 75 | 500 |
| ## | GAGAGC | 50 | 4.000 | 25 | 10 | 25 | 50 | 500 |
| ## | GAGTCT | 25 | 2.000 | 25 | 10 | 25 | 75 | 500 |
| ## | GAGTGA | 25 | 4.000 | 25 | 10 | 25 | 75 | 500 |
| ## | ACACAG | 100 | 0.000 | 0 | 10 | 25 | 25 | 500 |
| ## | ACATCT | 50 | 0.000 | 0 | 10 | 25 | 75 | 500 |
| ## | AGTCTC | 25 | 0.000 | 0 | 10 | 25 | 100 | 500 |
| ## | ATCGAT | 100 | 0.000 | 0 | 10 | 25 | 25 | 500 |
| ## | CACAGC | 50 | 0.000 | 0 | 10 | 25 | 75 | 500 |
| ## | CACTGA | 25 | 0.000 | 0 | 10 | 25 | 100 | 500 |
| ## | CGATAC | 100 | 0.000 | 0 | 10 | 25 | 25 | 500 |
| ## | CTGCGT | 50 | 0.000 | 0 | 10 | 25 | 75 | 500 |
| ## | GAGATA | 25 | 0.000 | 0 | 10 | 25 | 100 | 500 |
| ## | GCTAGC | 100 | 0.125 | 25 | 0 | 0 | 25 | 500 |
| ## | GCTATA | 100 | 0.250 | 25 | 0 | 0 | 25 | 500 |
| ## | GCTCAG | 100 | 0.500 | 25 | 0 | 0 | 25 | 500 |
| ## | GCTCGT | 100 | 1.000 | 25 | 0 | 0 | 25 | 500 |
| ## | GCTCTC | 100 | 2.000 | 25 | 0 | 0 | 25 | 500 |
| ## | GCTGAT | 100 | 4.000 | 25 | 0 | 0 | 25 | 500 |
| ## | GCTACG | 100 | 0.000 | 0 | 0 | 0 | 50 | 500 |
| ## | PRIMERS_RATIO | | | | | | | |
| ## | ACACGT | 640 | | | | | | |
| ## | ACACTC | 320 | | | | | | |
| ## | ACATGA | 320 | | | | | | |
| ## | ACAGAT | 160 | | | | | | |
| ## | AGTACG | 160 | | | | | | |
| ## | AGTGAT | 160 | | | | | | |
| ## | ACAGCA | 80 | | | | | | |
| ## | AGTAGC | 80 | | | | | | |
| ## | AGTGCA | 80 | | | | | | |
| ## | ATCGCA | 80 | | | | | | |
| ## | ACAGTG | 40 | | | | | | |
| ## | AGTATA | 40 | | | | | | |
| ## | AGTGTG | 40 | | | | | | |
| ## | ATCGTG | 40 | | | | | | |
| ## | CACATA | 40 | | | | | | |
| ## | ACATAC | 20 | | | | | | |
| ## | AGTCAG | 20 | | | | | | |
| ## | ATCACG | 20 | | | | | | |

| | |
|-------------------------|----------|
| ## ATCTAC | 20 |
| ## CACGAT | 20 |
| ## CGACAG | 20 |
| ## AGTCGT | 10 |
| ## ATCAGC | 10 |
| ## ATCTCT | 10 |
| ## CACGCA | 10 |
| ## CGACGT | 10 |
| ## CGATCT | 10 |
| ## ATCATA | 5 |
| ## ATCTGA | 5 |
| ## CACGTG | 5 |
| ## CGACTC | 5 |
| ## CGATGA | 5 |
| ## CTGCTC | 5 |
| ## CACACG | 2.5 |
| ## CACTAC | 2.5 |
| ## CGAGAT | 2.5 |
| ## CTGACG | 2.5 |
| ## CTGTAC | 2.5 |
| ## GAGCAG | 2.5 |
| ## CACTCT | 1.25 |
| ## CGAGCA | 1.25 |
| ## CTGAGC | 1.25 |
| ## CTGTCT | 1.25 |
| ## GAGCGT | 1.25 |
| ## CGAGTG | 0.625 |
| ## CTGATA | 0.625 |
| ## CTGTGA | 0.625 |
| ## GAGCTC | 0.625 |
| ## CTGCAG | 0.3125 |
| ## GAGACG | 0.3125 |
| ## GAGTAC | 0.3125 |
| ## GAGAGC | 0.15625 |
| ## GAGTCT | 0.15625 |
| ## GAGTGA | 0.078125 |
| ## ACACAG no_RT_PRIMERS | |
| ## ACATCT no_RT_PRIMERS | |
| ## AGTCTC no_RT_PRIMERS | |
| ## ATCGAT no_RT_PRIMERS | |
| ## CACAGC no_RT_PRIMERS | |
| ## CACTGA no_RT_PRIMERS | |
| ## CGATAC no_RT_PRIMERS | |
| ## CTGCGT no_RT_PRIMERS | |
| ## GAGATA no_RT_PRIMERS | |
| ## GCTAGC | 80 |
| ## GCTATA | 40 |
| ## GCTCAG | 20 |
| ## GCTCGT | 10 |
| ## GCTCTC | 5 |
| ## GCTGAT | 2.5 |
| ## GCTACG no_RT_PRIMERS | |

Create a CAGEexp object and load expression data

Number of sequencing reads extracted per sample

```
data.frame(colData(myCAGEexp)[,"extracted",drop=F])
```

```
##      extracted
## ACACGT      39556
## ACACTC     159844
## ACATGA      32474
## ACAGAT     138710
## AGTACG      14180
## AGTGAT     168934
## ACAGCA     158662
## AGTAGC      56598
## AGTGCA     428768
## ATCGCA      43925
## ACAGTG     218169
## AGTATA      15195
## AGTGTG      16988
## ATCGTG       5852
## CACATA      49911
## ACATAC     252993
## AGTCAG      18231
## ATCACG      11697
## ATCTAC      61931
## CACGAT      14078
## CGACAG       5787
## AGTCGT      42894
## ATCAGC    1053249
## ATCTCT       7558
## CAGGCA     168925
## CGACGT      50902
## CGATCT      26817
## ATCATA      17945
## ATCTGA     198657
## CACGTG      30261
## CGACTC     151581
## CGATGA       3492
## CTGCTC      12831
## CACACG    1333918
## CACTAC     551765
## CGAGAT     270601
## CTGACG      15719
## CTGTAC       7689
## GAGCAG      11251
## CACTCT     399697
## CGAGCA     135396
## CTGAGC      31942
## CTGTCT     231518
## GAGCGT      11654
## CGAGTG      43603
## CTGATA      45389
```

```
## CTGTGA      29711
## GAGCTC      19222
## CTGCAG      92305
## GAGACG      10894
## GAGTAC       917
## GAGAGC      92220
## GAGTCT       2006
## GAGTGA      17132
## ACACAG      42054
## ACATCT      33745
## AGTCTC       3059
## ATCGAT      10403
## CACAGC       8329
## CACTGA       6506
## CGATAC       1633
## CTGCGT       1318
## GAGATA       717
## GCTAGC       1895
## GCTATA       1206
## GCTCAG       2855
## GCTCGT       3098
## GCTCTC       4055
## GCTGAT      23130
## GCTACG       915
```

CTSS ANALYSIS

Number of nanoCAGE tags mapping at CTSS positions in each group of samples

```
## Figures not displayed on the html/pdf output
plotReverseCumulatives(myCAGEexp[,1:7], onePlot = TRUE)

plotReverseCumulatives(myCAGEexp[,8:14], onePlot = TRUE)

plotReverseCumulatives(myCAGEexp[,15:21], onePlot = TRUE)

plotReverseCumulatives(myCAGEexp[,22:28], onePlot = TRUE)

plotReverseCumulatives(myCAGEexp[,29:35], onePlot = TRUE)

plotReverseCumulatives(myCAGEexp[,36:42], onePlot = TRUE)

#plotReverseCumulatives(myCAGEexp$RNA == "10", onePlot = TRUE)
#plotReverseCumulatives(myCAGEexp[,50:56], onePlot = TRUE)
#plotReverseCumulatives(myCAGEexp[,57:63], onePlot = TRUE)
#plotReverseCumulatives(myCAGEexp[,64:70], onePlot = TRUE)
```

Number of nanoCAGE tags mapping at CTSS positions in each sample

```
(myCAGEexp$l1 <- colSums(CTSStagCountDf(myCAGEexp) > 0))

## ACACGT AACTC ACATGA ACAGAT AGTACG AGTGAT ACAGCA AGTAGC AGTGCA ATCGCA
```

```
## 3608 12843 3441 10063 1498 14452 8617 4358 23072 3661
## ACAGTG AGTATA AGTGTG ATCGTG CACATA ACATAC AGTCAG ATCACG ATCTAC CACGAT
## 12286 1646 1759 633 4004 11906 1607 1244 4741 1196
## CGACAG AGTCGT ATCAGC ATCTCT CACGCA CGACGT CGATCT ATCATA ATCTGA CACGTG
## 453 2671 40564 711 6328 2572 1882 1352 11459 1469
## CGACTC CGATGA CTGCTC CACACG CACTAC CGAGAT CTGACG CTGTAC GAGCAG CACTCT
## 5258 266 984 31921 16080 9138 878 484 775 12501
## CGAGCA CTGAGC CTGTCT GAGCGT CGAGTG CTGATA CTGTGA GAGCTC CTGCAG GAGACG
## 4243 1526 4360 160 1566 2119 1637 1231 2991 650
## GAGTAC GAGAGC GAGTCT GAGTGA ACACAG ACATCT AGTCTC ATCGAT CACAGC CACTGA
## 96 3592 182 960 4974 4174 525 1675 934 1015
## CGATAC CTGCGT GAGATA GCTAGC GCTATA GCTCAG GCTCGT GCTCTC GCTGAT GCTACG
## 145 148 95 345 167 422 410 195 1309 218
```

Create CTSS clusters

```
#clusterCTSS(myCAGEexp, thresholdIsTpm = FALSE, useMulticore = TRUE, nrPassThreshold = 2, removeSinglet
#cumulativeCTSSdistribution(myCAGEexp, clusters = "tagClusters")
##, use multicore = TRUE)
#quantilePositions(myCAGEexp, clusters = "tagClusters", qLow = 0.1, qUp = 0.9, useMulticore = TRUE)
```

Annotation with GENCODE

Collect Gencode annotations and gene symbols via AnnotationHub.

```
ah <- AnnotationHub::AnnotationHub()
ah["AH49547"]

## AnnotationHub with 1 record
## # snapshotDate(): 2017-04-25
## # names(): AH49547
## # $dataprovder: Gencode
## # $species: Mus musculus
## # $rdataclass: GRanges
## # $rdatadateadded: 2015-08-14
## # $title: gencode.vM6.basic.annotation.gff3.gz
## # $description: Gene annotations on reference chromosomes from Gencode
## # $taxonomyid: 1090
## # $genome: GRCm38
## # $sourcetype: GFF
## # $sourceurl: ftp://ftp.sanger.ac.uk/pub/gencode/Gencode_mouse/release_...
## # $sourcesize: 20384812
## # $tags: c("gencode", "vM6", "basic", "annotation", "gff3")
## # retrieve record with 'object[["AH49547"]]'
```

Annotate the genomic ranges of the tagCountMatrix SummarizedExperiment.

```
annotateCTSS(myCAGEexp, ah[["AH49547"]])
#annotateConsensusClusters(myCAGEexp, ah[["AH49547"]])
#consensusClustersSE(myCAGEexp)
#consensusClustersGR(myCAGEexp)
```

Make a gene expression table (not really required now).

```
CTSSStoGenes(myCAGEexp)
#CTSScoordinatesGR(myCAGEexp)
```

Save myCAGEexp file.


```
saveRDS(myCAGEexp, "myCAGEexp.Rds")
```

QC PLOTS

Boxplots


Extracted reads

```
par(mar=c(5,5,2,2), cex.main = 1, font.main = 2)
boxplot(extracted ~ barcode_ID, xlab = "Barcode ID", ylab = "Extracted reads", data = colData(myCAGEexp))
```




```
par(mar=c(5,5,2,2), cex.main = 1, font.main = 2)
boxplot(extracted ~ RT_PRIMERS, ylab = "Extracted reads", xlab = "RT_PRIMERS concentration (uM)", data = colData(myCAGEexp))
```

Number of extracted reads for different RT_PRIMERS concentrations




```
par(mar=c(5,5,2,2), cex.main = 1, font.main = 2)
boxplot(extracted ~ TSO, ylab = "Extracted reads", xlab = "TSO concentration (uM)", data = colData(myCA
```

Number of extracted reads for different TSO concentrations




```
par(mar=c(7,5,2,2), cex.main = 1, font.main = 2)
boxplot(extracted ~ PRIMERS_RATIO, ylab = "Extracted reads", xlab = "Primers ratio", data = colData(myC
```




Artefacts


```
par(mar=c(5,5,2,2), cex.main = 1, font.main = 2)
boxplot(tagdust / extracted ~ barcode_ID, xlab = "Barcode ID", ylab = "Artefacts/extracted reads", data
```


```
par(mar=c(5,5,2,2), cex.main = 1, font.main = 2)
boxplot(tagdust / extracted ~ RT_PRIMERS, ylab = "Artefacts/extracted reads", xlab = "RT_PRIMERS concentration (uM)", main = "Ratio artefacts/extracted for different RT_PRIMERS concentrations")
```



```
par(mar=c(5,5,2,2), cex.main = 1, font.main = 2)
boxplot(tagdust / extracted ~ TSO, ylab = "Artefacts/extracted reads", xlab = "TSO concentration (uM)",
```




```
par(mar=c(7,5,2,2), cex.main = 1, font.main = 2)
boxplot(tagdust / extracted ~ PRIMERS_RATIO, ylab = "Artefacts/extracted reads", xlab = "Primers ratio"
```




Counts


```
par(mar=c(5,5,2,2), cex.main = 1, font.main = 2)
boxplot(counts ~ barcode_ID, xlab = "Barcode ID", ylab = "Counts", data = colData(myCAGEexp), cex.axis = 1.5)
```




```
par(mar=c(5,5,2,2), cex.main = 1, font.main = 2)
boxplot(counts ~ RT_PRIMERS, ylab = "Counts", xlab = "RT_PRIMERS concentration (uM)", data = colData(my
```



```
par(mar=c(5,5,2,2), cex.main = 1, font.main = 2)
boxplot(counts ~ TSO, ylab = "Counts", xlab = "TSO concentration (uM)", data = colData(myCAGEexp), cex.
```




```
par(mar=c(7,5,2,2), cex.main = 1, font.main = 2)
boxplot(counts ~ PRIMERS_RATIO, ylab = "Counts", xlab = "Primers ratio", data = colData(myCAGEexp), cex.
```




Genes


```
par(mar=c(5,5,2,2), cex.main = 1, font.main = 2)
boxplot(genes ~ barcode_ID, xlab = "Barcode ID", ylab = "Genes", data = colData(myCAGEexp), cex.axis = 0.8)
```




```
par(mar=c(5,5,2,2), cex.main = 1, font.main = 2)
boxplot(genes ~ RT_PRIMERS, ylab = "Genes", xlab = "RT_PRIMERS concentration (uM)", data = colData(myCA
```



```
par(mar=c(5,5,2,2), cex.main = 1, font.main = 2)
boxplot(genes ~ TSO, ylab = "Genes", xlab = "TSO concentration (uM)", data = colData(myCAGEexp), cex.a
```




```
par(mar=c(7,5,2,2), cex.main = 1, font.main = 2)
boxplot(genes ~ PRIMERS_RATIO, ylab = "Genes", xlab = "Primers ratio", data = colData(myCAGEexp), cex.a
```

Dotplots


Extracted reads

```
dotsize <- 25000
ggplot(colData(myCAGEexp) %>% data.frame, aes(x=PRIMERS_RATIO, y=extracted)) +
  stat_summary(fun.y=mean, fun.ymin=mean, fun.ymax=mean, geom="crossbar", color="gray") +
  geom_dotplot(aes(fill=group), binaxis='y', binwidth=1.5, dotsize=dotsize, stackdir='center') + theme_l
  xlab("Primers ratio\n (TSO/RT_primers)") +
  labs(title = "Number of sequencing reads extracted") +
  coord_flip()
```



Artefacts


```
dotsize <- 5000
ggplot(colData(myCAGEexp) %>% data.frame, aes(x=PRIMERS_RATIO, y=tagdust)) +
  stat_summary(fun.y=mean, fun.ymin=mean, fun.ymax=mean, geom="crossbar", color="gray") +
  geom_dotplot(aes(fill=group), binaxis='y', binwidth=1, dotsize=dotsize, stackdir='center') + theme_bw()
xlab("Primers ratio\n (TSO/RT_primers)") +
labs(title = "Artefacts reads per sample") +
coord_flip()
```



Counts


Gene counts

```
dotsize <- 100
ggplot(colData(myCAGEexp) %>% data.frame, aes(x=PRIMERS_RATIO, y=genes)) +
  stat_summary(fun.y=mean, fun.ymin=mean, fun.ymax=mean, geom="crossbar", color="gray") +
  geom_dotplot(aes(fill=group), binaxis='y', binwidth=1, dotsize=dotsize, stackdir='center') + theme_bw()
xlab("Primers ratio\n (TSO/RT_primers)") +
labs(title = "Average number of genes detected in each group") +
coord_flip()
```



Transcript counts


```
dotsize <- 4000
ggplot(colData(myCAGEexp) %>% data.frame, aes(x=PRIMERS_RATIO, y=counts)) +
  stat_summary(fun.y=mean, fun.ymin=mean, fun.ymax=mean, geom="crossbar", color="gray") +
  geom_dotplot(aes(fill=group), binaxis='y', binwidth=1, dotsize=dotsize, stackdir='center') + theme_bw()
xlab("Primers ratio\n (TSO/RT_primers)") +
  labs(title = "Average number of molecules detected in each group") +
  coord_flip()
```



Barplots


Quality control by sample

```
plotAnnot(myCAGEexp, "qc", "Quality control by sample", "barcode_ID", normalise = F)
```



Annotations by sample

```
plotAnnot(myCAGEexp, "counts", "Annotations by sample", "barcode_ID", normalise = F)
```




Rarefaction

```
rar1 <- hanabi(CTSStagCountDF(myCAGEexp), from = 0)
rarc <- hanabi(assay(consensusClustersSE(myCAGEexp)) %>% as.data.frame, from = 0)
rarg <- hanabi(assay(GeneExpSE(myCAGEexp)) %>% as.data.frame, from = 0)
save(rar1, rarg, file="rar.Rda")
```

Plot TSS discovery

```
hanabiPlot(rar1, ylab='number of TSS detected', xlab='number of unique molecule counts', main=paste("TSS"))
```

TSS discovery




Plot Cluster discovery

```
#hanabiPlot(rarc, ylab='number of CTSS clusters detected', xlab='number of unique molecule counts', main='Cluster discovery')
```

Plot Gene discovery

```
hanabiPlot(rarg, ylab='number of genes detected', xlab='number of unique molecule counts', main=paste("Gene discovery", rarg))
```



Correlation of CTSS clusters expression across samples (sample distance)

Create the correlation matrix

```
#c <- assay(consensusClustersSE(myCAGEexp)) %>% as.data.frame
#cor_clusters <- cor(log1p(c))
#cor_clusters[cor_clusters == 1] <- NA
```

Heatmap

```
#col <- colorRampPalette( rev(brewer.pal(9, "RdBu")) )(255)
#ann_col <- data.frame(row.names = rownames(colData(myCAGEexp)), protocol = myCAGEexp$protocol)
#ann_col$protocol <- sub("CONTROL", "CTL", ann_col$protocol)
#ann_row <- data.frame(row.names = rownames(colData(myCAGEexp)), step = myCAGEexp$step)
#ann_row$step <- sub("0", "CTL", ann_row$step)
#pheatmap::pheatmap(cor_clusters, show_colnames = F, annotation_row = ann_row, annotation_col = ann_col)
```


Principal Component Analysis (PCA)

Define PCA axis based on correlation matrix

```
#cor_clusters <- cor(log1p(c))  
#PCA <- prcomp(cor_clusters, scale. = TRUE)  
#dfdf <- stats::summary.prcomp(PCA)$importance[2, ]
```

List of principal components identified ranked by % of explained variability:

PCA Plot 1: PC1 vs. PC2 annotated by PRIMERS_RATIO

```
#ggbiplot::ggbiplot(PCA, choices=c(1,2), obs.scale = 1, var.scale = 1, groups = colData(myCAGEexp)$resu
```

Plate maps

Create plate object

```
plate <- as.data.frame(colData(myCAGEexp))
```

```
plateMap <- function(x, title) {  
  platetools::raw_map(plate[[x]], well=plate$Well, plate="384") +  
  ggtitle(title) +  
  viridis::scale_fill_viridis(breaks = unique(plate[[x]]))  
}  
  
plateMapLog <- function(x, title) {  
  platetools::raw_map(plate[[x]], well=plate$Well, plate="384") +  
  ggtitle(title) +  
  viridis::scale_fill_viridis(breaks = unique(plate[[x]]), trans = "log")  
}
```

RT primers


```
(plot_RT <- plateMapLog("RT_PRIMERS", "RT primer concentration"))
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.  
## Are you sure argument 'plate' is correct for the number of wells in your data?  
## note: Default is a 96-well plate.
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.  
## Are you sure argument 'plate' is correct for the number of wells in your data?  
## note: Default is a 96-well plate.
```

```
## Warning: Transformation introduced infinite values in discrete y-axis
```

RT primer concentration




TSO

```
(plot_TSO <- plateMapLog("TSO", "TSO concentration"))
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

TSO concentration




Ratio TSO / RT

```
(plot_TSO_RT_ratio <- platetools::raw_map(plate$TSO / plate$RT_PRIMERS, well=plate$Well, plate="384") +
  ggtitle("TSO / RT primer concentration") +
  viridis::scale_fill_viridis(breaks = unique(plate$TSO / plate$RT_PRIMERS), trans = "log"))
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

TSO / RT primer concentration




Extracted

```
(plot_extracted <- plateMapLog("extracted", "Extracted reads"))
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

Extracted reads




Artefacts

```
(plot_artefacts <- plateMapLog("tagdust", "Artefacts"))
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

Artefacts




rDNA

```
(plot_rDNA <- plateMapLog("rdna", "Reads mapping to rDNA"))
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

Reads mapping to rDNA




Genes

```
(plot_genes <- plateMapLog("genes", "Genes detected"))
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

Genes detected




Counts

```
(plot_counts <- plateMapLog("counts", "Molecules detected"))
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```


Molecules detected




Promoters

```
(plot_promoters <- plateMapLog("promoter", "Promoters detected"))
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

```
## Warning: Invalid plate selection. The data given has more rows then number of wells.
## Are you sure argument 'plate' is correct for the number of wells in your data?
## note: Default is a 96-well plate.
```

Promoters detected




Grand summary plots

Extracted. artefacts, rDNA

```
ggpubr::ggarrange(ncol = 3, nrow = 2, plot_RT, plot_TSO, plot_TSO_RT_ratio, plot_extracted, plot_artefacts)
```


Warning: Transformation introduced infinite values in discrete y-axis



Molecules, genes and promoters

```
ggpubr::ggarrange(ncol = 3, nrow = 2, plot_RT, plot_TSO, plot_TSO_RT_ratio, plot_counts, plot_genes, plot_promoters)
```

```
## Warning: Transformation introduced infinite values in discrete y-axis
```



```
sessionInfo()
```

```
## R version 3.4.0 (2017-04-21)
## Platform: x86_64-pc-linux-gnu (64-bit)
## Running under: Debian GNU/Linux 9 (stretch)
##
## Matrix products: default
## BLAS: /usr/lib/libblas/libblas.so.3.7.0
## LAPACK: /usr/lib/lapack/liblapack.so.3.7.0
##
## locale:
##  [1] LC_CTYPE=C.UTF-8      LC_NUMERIC=C           LC_TIME=C.UTF-8
##  [4] LC_COLLATE=C.UTF-8    LC_MONETARY=C.UTF-8    LC_MESSAGES=C.UTF-8
##  [7] LC_PAPER=C.UTF-8      LC_NAME=C              LC_ADDRESS=C
## [10] LC_TELEPHONE=C        LC_MEASUREMENT=C.UTF-8 LC_IDENTIFICATION=C
##
## attached base packages:
## [1] stats4      parallel  stats      graphics  grDevices  utils      datasets
## [8] methods    base
##
## other attached packages:
##  [1] bindrcpp_0.2                vegan_2.4-4
##  [3] lattice_0.20-35             permute_0.9-4
##  [5] reshape_0.8.7               SummarizedExperiment_1.6.5
##  [7] DelayedArray_0.2.7          matrixStats_0.52.2
##  [9] Biobase_2.36.2              MultiAssayExperiment_1.2.1
## [11] plyr_1.8.4                  magrittr_1.5
## [13] RColorBrewer_1.1-2          plots_3.0.1
## [15] ggplot2_2.2.1               data.table_1.10.4-1
## [17] CAGEr_1.21.0                BSgenome.Mmusculus.UCSC.mm9_1.4.0
## [19] BSgenome_1.44.2             rtracklayer_1.36.4
## [21] Biostings_2.44.2            XVector_0.16.0
## [23] GenomicRanges_1.28.4        GenomeInfoDb_1.12.3
## [25] IRanges_2.10.3              S4Vectors_0.14.5
## [27] BiocGenerics_0.22.1
##
## loaded via a namespace (and not attached):
##  [1] nlme_3.1-131                bitops_1.0-6
```

| | |
|------------------------------|-------------------------------|
| ## [3] bit64_0.9-7 | httr_1.3.1 |
| ## [5] rprojroot_1.3-2 | UpSetR_1.3.3 |
| ## [7] tools_3.4.0 | backports_1.1.2 |
| ## [9] platetools_0.0.2 | R6_2.2.2 |
| ## [11] KernSmooth_2.23-15 | DBI_0.7 |
| ## [13] lazyeval_0.2.1 | mgcv_1.8-17 |
| ## [15] colorspace_1.3-2 | gridExtra_2.3 |
| ## [17] curl_3.1 | bit_1.1-12 |
| ## [19] compiler_3.4.0 | VennDiagram_1.6.18 |
| ## [21] labeling_0.3 | caTools_1.17.1 |
| ## [23] scales_0.5.0 | stringr_1.2.0 |
| ## [25] digest_0.6.13 | Rsamtools_1.28.0 |
| ## [27] rmarkdown_1.8 | pkgconfig_2.0.1 |
| ## [29] htmltools_0.3.6 | rlang_0.1.2 |
| ## [31] RSQLite_2.0 | VGAM_1.0-4 |
| ## [33] BiocInstaller_1.26.1 | shiny_1.0.5 |
| ## [35] bindr_0.1 | BiocParallel_1.10.1 |
| ## [37] gtools_3.5.0 | dplyr_0.7.4 |
| ## [39] RCurl_1.95-4.9 | GenomeInfoDbData_0.99.0 |
| ## [41] futile.logger_1.4.3 | smallCAGEqc_0.12.2.99999 |
| ## [43] Matrix_1.2-10 | Rcpp_0.12.12 |
| ## [45] munsell_0.4.3 | viridis_0.4.0 |
| ## [47] stringi_1.1.6 | yaml_2.1.16 |
| ## [49] MASS_7.3-47 | zlibbioc_1.22.0 |
| ## [51] AnnotationHub_2.8.3 | blob_1.1.0 |
| ## [53] grid_3.4.0 | gdata_2.18.0 |
| ## [55] shinydashboard_0.6.1 | cowplot_0.9.2 |
| ## [57] splines_3.4.0 | knitr_1.18 |
| ## [59] beanplot_1.2 | ggpubr_0.1.6 |
| ## [61] reshape2_1.4.3 | codetools_0.2-15 |
| ## [63] futile.options_1.0.0 | XML_3.98-1.9 |
| ## [65] glue_1.2.0 | evaluate_0.10.1 |
| ## [67] lambda.r_1.2 | httpuv_1.3.5 |
| ## [69] gtable_0.2.0 | purrr_0.2.4 |
| ## [71] tidyr_0.7.1 | assertthat_0.2.0 |
| ## [73] mime_0.5 | xtable_1.8-2 |
| ## [75] viridisLite_0.2.0 | tibble_1.3.4 |
| ## [77] som_0.3-5.1 | AnnotationDbi_1.38.2 |
| ## [79] memoise_1.1.0 | GenomicAlignments_1.12.2 |
| ## [81] cluster_2.0.6 | interactiveDisplayBase_1.14.0 |