

# **Labcyte 525 demo test results**

Charles Plessy, RIKEN, June 14<sup>th</sup>, 2018

# Parameter space exploration of a biochemical reaction

Using the the nanoliter-scale liquid handling platform “Labcyte Echo” platform, we assembled more than 3,000 reverse-transcription (RT) reactions in eight 384-well plates. We then assessed the RT efficiency by sequencing them with the nanoCAGE method. The RT reactions were made in quadruplicates and differed in the combination of the following reaction parameters:


- Template switching oligonucleotide (TSO): 9 different concentrations ranging from 0.6 to 160  $\mu\text{M}$ .
- Reverse-transcription primer (RTP): 6 different concentration ranging from 1 to 24  $\mu\text{M}$ .
  - RNA quantity: 6 different starting amounts from 1 pg to 100 ng.
  - Enzyme type (SuperScript III or IV)

This makes 2,592 reactions, to which we added 480 negative controls.

# Rationale

- In some projects, we need to know the range of reagent concentrations in which the reaction works, because, because of experimental constraints limiting them. (For instance highest concentrations may be too viscous to handle, or lowest concentrations may be hard to produce because of adsorption losses, etc.)
- Some of our past benchmarks are hard to replicate accurately across long periods of time (multiple years), when material, operators, reagent batches and laboratory location have fully turned over. A standardised automation platform may be the key for long-term replication.
- With a systematic approach, we hope to comprehensively characterise a biochemical system to the point that iterative approaches over years are not necessary anymore. (Without the high resolution reached in this work, developments are made by small incremental changes scattered over years – also known as “protocol drift”).


# Computer-assisted experiment design: example reagent maps from the 384-well plates that we prepared



*Above: different randomisation on different replicates*


Computer-generated randomisation plays an important role in our study: to eliminate positional biases, the well coordinates of reactions with the same reagent parameters are randomly chosen. Moreover, the barcode identifiers (for multiplexing next-generation sequencing) are also randomised. Lastly, the position of the negative controls (for instance the grey squares in the “RNA mass” panel) are randomised, thus generating a plate fingerprint that prevents accidental sample swaps. Here, the RNA mass parameter was not randomised because different masses needed different numbers of PCR cycles in the nanoCAGE reactions. For practical purposes (multichannel pipetting and limited number of barcodes), each plate was divided in 6 interleaved “sextants”, but development of 384 different barcodes would entirely be possible.

# How to read the plots



# Too much but also too few oligonucleotides cause high amounts of artefacts


Amount of oligonucleotide artefacts (% of extracted reads)



# Identification of non-linear patterns

(when looking at the amount of ribosomal RNA detected – lower is better)

Fraction of reads aligning to rRNA sequences (% of non-tagdust extracted reads)




At lower resolution (less points on a narrower dynamic range), we could have concluded that there is a monotonic response followed by a plateau.

# Counter-intuitive findings: more TSO artefacts at low TSO concentrations !

Strand Invasion (% of molecule counts)

One outlier removed (1 pg RNA, SSIII, 1  $\mu$ M RTP, > 30% strand invasion).



With such experiments, we gain a mechanistic insight.

Linking promoters to functional transcripts in small samples with nanoCAGE and CAGEscan

Charles Plessy<sup>1,10</sup>, Nicolas Bertin<sup>1,10</sup>, Hazuki Takahashi<sup>1,10</sup>, Roberto Simone<sup>2,10</sup>, Md. Salimullah<sup>1</sup>, Timo Lassmann<sup>1</sup>, Morana Vitezic<sup>1,3</sup>, Jessica Severin<sup>1</sup>, Signe Olivarius<sup>1,8</sup>, Dejan Lazarevic<sup>2</sup>, Nadine Hornig<sup>1</sup>, Valerio Orlando<sup>2</sup>, Ian Bell<sup>1</sup>, Hui Gao<sup>1</sup>, Jacqueline Dumais<sup>4</sup>, Philipp Kapranov<sup>5,9</sup>, Huihui Wang<sup>6</sup>, Carrie A. Davis<sup>6</sup>, Thomas R. Gingeras<sup>6</sup>, Jun Kawai<sup>1</sup>, Carsten O. Daub<sup>1</sup>, Yoshihide Hayashizaki<sup>1</sup>, Stefano Gustincich<sup>2</sup>, and Piero Carninci<sup>1</sup>

<sup>1</sup>RIKEN Yokohama Institute, Omics Science Center (OSC), Yokohama, Japan

<sup>2</sup>Sector of Neurobiology, Scuola Internazionale Superiore di Studi Avanzati (SISSA), Trieste, Italy


<sup>3</sup>Department of Cell and Molecular Biology (CMB), Karolinska Institutet, Stockholm, Sweden

<sup>4</sup>Affymetrix Inc., Santa Clara, California, USA

<sup>5</sup>Cold Spring Harbor Laboratory, Genome Center, Woodbury, New York, USA

<sup>6</sup>Watson School of Biological Sciences, Cold Spring Harbor Laboratory, Cold Spring Harbor, New York, USA

<sup>7</sup>Dulbecco Telethon Institute, Istituto Di Ricovero e Cura a Carattere Scientifico (IRCCS) Fondazione Santa Lucia, Rome, Italy



Published online 24 November 2012

Nucleic Acids Research, 2013, Vol. 41, No. 3 e64  
doi:10.1093/nar/gks1128

Suppression of artifacts and barcode bias in high-throughput transcriptome analyses utilizing template switching

Dave T. P. Tang<sup>1</sup>, Charles Plessy<sup>1</sup>, Md Salimullah<sup>1</sup>, Ana Maria Suzuki<sup>1</sup>, Raffaella Calligaris<sup>2</sup>, Stefano Gustincich<sup>2</sup> and Piero Carninci<sup>1\*</sup>

<sup>1</sup>Omics Science Center, RIKEN Yokohama Institute, 1-7-22 Suehiro-cho, Tsurumi-ku, Yokohama, Kanagawa 230-0045, Japan and <sup>2</sup>Sector of Neurobiology, International School for Advanced Studies (SISSA), via Bonomea 265, 34134 Trieste, Italy

Received March 13, 2012; Revised September 27, 2012; Accepted October 23, 2012





# High TS oligos is almost always good

(thus, we should not fear that beads-based method would introduce too much of them)

## Promoter rate

(% of molecule counts after SI removal)



← alternative representation

# Conclusions and perspectives

- Further directions
  - Put in perspective with software parameter optimisation and AI.
  - Try other graphical representations ? (3D plots, ratios of parameters, ...)
  - Explore more parameters (dNTP and divalent concentrations)
  - Further automate (3D printed sample collector, library preparation on Labcyte platform, etc.)
- Conclusions regarding reverse-transcription
  - In general, high [TSO] and low [RTP] is better.
  - But [RTP] should be matched to [TSO].
  - Optima are different at single-cell scale (10 pg) and bulk scale (100 ng).
- Perspectives
  - Optimise other reactions (ligations, template-switching of miRNAs, downstream PCR steps, ...)
  - Further automation of experiment design.

# Acknowledgements

- RIKEN:
  - Sachi Kato
  - Stéphane Poulain
  - Aki Minoda
  - Piero Carninci
- Labcyte
  - Iris Chen
  - Hiro Ishida
- Kiko-tech
  - Hisakazu Komina